Experiment Name: Determination of Solids in Water as per Indian Standards

Objective:

To measure the various types of solids (total solids, suspended solids, dissolved solids, and volatile solids) in water as per IS 3025 (Part 15): 1984.

Apparatus and Materials

1. Glassware:

- o Evaporating dishes (porcelain or glass).
- o Gooch crucibles or Whatman filter papers (for suspended solids).
- o Beakers (250 mL).
- 2. **Drying Oven**: Capable of maintaining a temperature of 103–105°C.
- 3. **Muffle Furnace**: Capable of reaching 550°C for volatile solids determination.
- 4. **Analytical Balance**: Sensitivity of 0.1 mg.
- 5. Vacuum Pump or Filtration Apparatus: For filtration of suspended solids.
- 6. **Desiccator**: To cool samples to room temperature in a moisture-free environment.
- 7. **Distilled Water**: For rinsing and sample preparation.
- 8. Glass Fiber Filter Discs (or equivalent): For separation of suspended solids.

Types of Solids

- 1. Total Solids (TS): Includes both dissolved and suspended solids.
- 2. Suspended Solids (SS): Solids retained on a filter after sample filtration.
- 3. **Dissolved Solids (DS)**: Solids that pass through the filter.
- 4. **Volatile Solids (VS)**: Portion of solids lost upon ignition at 550°C, indicative of organic content.

Procedure

A. Total Solids (TS)

1. **Preparation**:

- o Clean and dry an evaporating dish in a drying oven at 103–105°C.
- o Cool the dish in a desiccator and weigh it. Record the weight as W1.

2. Sample Measurement:

o Pour a measured volume (50–100 mL) of well-mixed water into the dish.

3. Evaporation and Drying:

- o Place the dish in the drying oven at 103–105°C for 24 hours.
- o Remove, cool in a desiccator, and weigh. Record the final weight as W2.
- 4. Calculation: Total Solids (mg/L)=(W2-W1)×1000/Volume of sample (mL)

B. Suspended Solids (SS)

1. **Preparation**:

o Dry a filter paper or glass fiber filter disc at 103–105°C. Cool and weigh it (W1).

2. Filtration:

- o Filter a measured volume (50–100 mL) of the water sample through the filter.
- o Rinse the filter with distilled water to remove dissolved solids.

3. Drying and Weighing:

- o Dry the filter at 103–105°C until a constant weight is obtained.
- o Cool in a desiccator and weigh (W2).
- **4.** Calculation: Suspended Solids (mg/L)=(W2-W1)×1000/Volume of sample (mL)

C. Dissolved Solids (DS)

1. Procedure:

 The dissolved solids are calculated as the difference between the total solids and suspended solids:

Dissolved Solids (mg/L) = Total Solids—Suspended Solids

D. Volatile Solids (VS)

1. **Preparation**:

- o Place the dried residue (from TS or SS determination) in a muffle furnace at 550°C for 15–30 minutes.
- o Allow the sample to cool in a desiccator and weigh it (W3).
- 2. Calculation: Volatile Solids (mg/L)=Total Solids-Ash (residue at 550°C)

Observation Table

Parameter	Sample Volume (mL)	Initial Weight (W ₁) (g)	Final Weight (W ₂ /W ₃) (g)	Weight Difference (g)	Concentration (mg/L)
Total Solids					
(TS)					
Suspended					
Solids (SS)					
Dissolved					
Solids (DS)					
Volatile Solids					
(VS)					

Precautions

- 1. Ensure consistent drying and cooling times for accurate results.
- 2. Avoid contamination of samples during handling and weighing.
- 3. Use freshly prepared or well-preserved samples.
- 4. Handle hot dishes and filters with care.

Results

The concentrations of total solids, suspended solids, dissolved solids, and volatile solids in the water sample are expressed in mg/L.

References

- IS 3025 (Part 15): 1984 Methods of Sampling and Test (Physical and Chemical) for Water and Wastewater: Determination of Solids.
- Relevant guidelines from the Central Pollution Control Board (CPCB), India.